Intended use

Orion Diagnostica Ultrasensitive CRP Control is intended for quality control of CRP assays measuring low CRP concentrations by a clinical chemistry analyser, QuikRead® 101 Instrument or QuikRead go® Instrument.
For *in vitro* diagnostic use.

Principles of the procedure

Use the control in the same way as unknown specimen would be used, following the instructions for the assay procedure.

Reagents

The control is human-based material. It should be handled in the same way as patient sample.

Preparation

The control is ready-to-use. Allow the material to stabilise for 30 minutes at room temperature (18...25°C). Mix well. Do not shake.

Storage

Stable until the expiry date when stored at 2...8°C. Opened vial: Stable for 2 months at 2...8°C. Do not use any vial that has visible evidence of microbial growth.

Expected values

The individual result should fall within the indicated range. Each laboratory should establish its own mean and precision data.

Ultrasensitive CRP Control, Cat. No. 68257	Deutsch
Ultrasensitive CRP Control, Kat. Nr. 68257	

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Warnings and precautions

- Reagents contain <0.1% sodium azide, which is not considered a harmful concentration. Sodium azide liberates toxic gas when in contact with acids. Azides may react with metal plumbing, forming explosive compounds. Azide build-up can be avoided by flushing with a large volume of water when disposing of the reagents.
- Do not use the product after the expiry date marked on the outer package.
- Do not exceed the stability periods for opened reagents.
- Ultrasensitive CRP Control contains material of human origin, which has been tested and found negative for human immunodeficiency virus antibodies (HIV 1 and 2), hepatitis C virus (HCV) and hepatitis B surface antigen (HbsAg). As a precaution, handle Ultrasensitive CRP Control as if capable of transmitting infectious agents.

Disposal

- Dispose of contents according to national and local law.
- Materials of the components
Paper: Instructions for use
Cardboard: Outer box
Glass: Control vials
Plastic: Vial caps
Rubber: Stopper
- When used in accordance with Good Laboratory Practice, good occupational hygiene and the instructions for use, the reagents supplied should not present a hazard to health.

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